# **Voll Clinic API**

This is a sample server to Voll Clinic. This is a doctors and patients CRUD and you can schedule some appointment with a doctor or specialty.

More information: <a href="https://helloreverb.com">https://helloreverb.com</a>
Contact Info: hello@helloreverb.com

Version: 1.0.0 Apache 2.0

http://springdoc.org

### Access

1.

### Methods

[ Jump to Models ]

**Table of Contents** 

### **AppointmentController**

- GET /api/appointment/{id}
- POST /api/appointment
- DELETE /api/appointment

### **DoctorController**

- POST /api/doctor
- DELETE /api/doctor/{id}
- GET /api/doctor
- GET /api/doctor/{id}
- POST /api/doctor/search
- PUT /api/doctor

### **PatientController**

- POST /api/patient
- DELETE /api/patient/{id}
- GET /api/patient
- GET /api/patient/{id}
- POST /api/patient/search
- PUT /api/patient

### **UserAuthentication**

- POST /api/auth/login
- POST /api/auth/register

# **AppointmentController**

GET /api/appointment/{id}

Get the appointment data with it's ID (getOne2)

### Path parameters

<u>Up</u>

id (required)
Path Parameter — format: int64

# Return type Appointment

### Example data

Content-Type: application/json

```
"doctor" : {
  "specialty": "ORTHOPEDY",
  "phone": "phone",
  "name" : "name",
  "id" : 1,
  "isActive" : true,
  "email" : "email",
  "dni" : "dni"
"date" : "2000-01-23T04:56:07.000+00:00",
"cancellationReason" : "PATIENT_CANCELED",
"patient" : {
  "address" : {
    "zipcode" : "zipcode",
    "address" : "address",
    "city" : "city",
    "county": "county",
    "complement" : "complement"
  "phone": "phone",
  "name" : "name",
  "id" : 6,
  "isActive" : true,
  "email" : "email",
  "dni" : "dni"
},
"id" : 0
```

### **Produces**

This API call produces the following media types according to the Accept request header; the media type will be conveyed by the Content-Type response header.

• \*/\*

# Responses

200

**OK** Appointment

# POST /api/appointment



Schedule one appointment with a doctor or a specialty (scheduleAppointment)

#### Consumes

This API call consumes the following media types via the Content-Type request header:

• application/json

### Request body

body <u>DataScheduleAppointment</u> (required) Body Parameter —

# Return type <u>DataAppointment</u>

### Example data

Content-Type: application/json

```
"doctor" : {
    "specialty": "ORTHOPEDY",
    "phone": "phone",
    "name": "name",
    "id" : 1,
    "isActive" : true,
    "email" : "email",
    "dni" : "dni"
  },
  "date" : "2000-01-23T04:56:07.000+00:00",
  "patient" : {
    "address" : {
      "zipcode" : "zipcode",
      "address" : "address",
      "city" : "city",
      "county": "county",
      "complement" : "complement"
    "phone": "phone",
    "name" : "name",
    "id" : 6,
    "isActive" : true,
    "email" : "email",
    "dni" : "dni"
  },
  "id" : 0
}
```

### **Produces**

This API call produces the following media types according to the Accept request header; the media type will be conveyed by the Content-Type response header.

• \*/\*

# Responses

200

**OK** <u>DataAppointment</u>

# DELETE /api/appointment



Cancel one appointment with a reason (scheduleAppointment1)

#### Consumes

This API call consumes the following media types via the Content-Type request header:

• application/json

### Request body

```
body <u>DataCancelAppointment</u> (required)
Body Parameter —
```

## Return type

String

### Example data

Content-Type: application/json

```
пп
```

### **Produces**

This API call produces the following media types according to the Accept request header; the media type will be conveyed by the Content-Type response header.

• \*/\*

# Responses

200

OK String

# **DoctorController**

POST /api/doctor

<u>Up</u>

(createDoctor1)

### Consumes

This API call consumes the following media types via the Content-Type request header:

• application/json

### Request body

```
body <u>DataDoctorRegister</u> (required)
Body Parameter —
```

# Return type

**DataDoctor** 

### Example data

Content-Type: application/json

```
"address" : {
    "zipcode" : "zipcode",
    "address" : "address",
    "city" : "city",
    "county" : "county",
    "complement" : "complement"
},
    "phone" : "phone",
    "name" : "name",
    "id" : 0,
    "email" : "email",
```

```
"dni" : "dni"
}
```

#### **Produces**

This API call produces the following media types according to the Accept request header; the media type will be conveyed by the Content-Type response header.

• \*/\*

# Responses 200

**OK DataDoctor** 

DELETE /api/doctor/{id}



(deactivateDoctor)

### Path parameters

id (required)

Path Parameter - format: int64

### Return type

String

### Example data

Content-Type: application/json

11 11

### **Produces**

This API call produces the following media types according to the Accept request header; the media type will be conveyed by the Content-Type response header.

• \*/\*

# Responses 200

200

**OK String** 

GET /api/doctor



(getAll1)

### Query parameters

pagination (required) Query Parameter —

# Return type

<u>PageDataDoctorList</u>

### Example data

Content-Type: application/json

```
{
    "number" : 5,
```

```
"size" : 1,
  "last" : true,
  "numberOfElements" : 2,
  "totalPages" : 0,
  "pageable" : {
    "paged" : true,
    "pageNumber" : 9,
    "offset" : 7,
    "pageSize" : 3,
    "unpaged" : true
  },
  "sort" : {
    "unsorted" : true,
    "sorted" : true,
    "empty" : true
  "content" : [ {
    "specialty": "ORTHOPEDY",
    "name" : "name",
    "id" : 5,
    "email" : "email",
    "dni" : "dni"
    "specialty" : "ORTHOPEDY",
    "name" : "name",
    "id" : 5,
    "email" : "email",
    "dni" : "dni"
  } ],
  "first" : true,
  "totalElements" : 6,
  "empty" : true
}
```

#### **Produces**

This API call produces the following media types according to the Accept request header; the media type will be conveyed by the Content-Type response header.

• \*/\*

# Responses 200

**OK** PageDataDoctorList

```
GET /api/doctor/{id}
(getOne1)
Path parameters
  id (required)
    Path Parameter — format: int64
```

### Return type <u>Doctor</u>

Example data

Content-Type: application/json

<u>Up</u>

```
{
  "specialty" : "ORTHOPEDY",
  "phone" : "phone",
  "name" : "name",
  "id" : 1,
  "isActive" : true,
  "email" : "email",
  "dni" : "dni"
}
```

#### **Produces**

This API call produces the following media types according to the Accept request header; the media type will be conveyed by the Content-Type response header.

• \*/\*

# Responses 200

**OK Doctor** 

# POST /api/doctor/search

<u>Up</u>

(searchByCriteria1)

### Consumes

This API call consumes the following media types via the Content-Type request header:

• application/json

## Request body

```
body <u>doctor_search_body</u> (required)
Body Parameter —
```

# Return type

### **PageDataDoctorList**

## Example data

Content-Type: application/json

```
"number" : 5,
"size" : 1,
"last" : true,
"numberOfElements" : 2,
"totalPages" : 0,
"pageable" : {
  "paged" : true,
  "pageNumber" : 9,
  "offset" : 7,
  "pageSize" : 3,
  "unpaged" : true
},
"sort" : {
  "unsorted" : true,
  "sorted" : true,
  "empty" : true
},
```

```
"content" : [ {
    "specialty" : "ORTHOPEDY",
    "name" : "name",
    "id" : 5,
    "email" : "email",
    "dni" : "dni"
}, {
    "specialty" : "ORTHOPEDY",
    "name" : "name",
    "id" : 5,
    "email" : "email",
    "dni" : "dni"
} ],
    "first" : true,
    "totalElements" : 6,
    "empty" : true
}
```

#### **Produces**

This API call produces the following media types according to the Accept request header; the media type will be conveyed by the Content-Type response header.

• \*/\*

# Responses

200

**OK** PageDataDoctorList

```
PUT /api/doctor
```

<u>Up</u>

(updateDoctor1)

### Consumes

This API call consumes the following media types via the Content-Type request header:

• application/json

### Request body

```
body <u>DataDoctorUpdate</u> (required)
Body Parameter —
```

### Return type

String

### Example data

Content-Type: application/json

```
пп
```

### **Produces**

This API call produces the following media types according to the Accept request header; the media type will be conveyed by the Content-Type response header.

• \*/\*

### Responses

200

OK String

# **PatientController**

POST /api/patient

<u>Up</u>

(createDoctor)

### Consumes

This API call consumes the following media types via the Content-Type request header:

• application/json

### Request body

```
body <u>DataPatientRegister</u> (required)
Body Parameter —
```

## Return type <u>DataPatient</u>

### Example data

Content-Type: application/json

```
{
   "address" : {
        "zipcode" : "zipcode",
        "address" : "address",
        "city" : "city",
        "county" : "county",
        "complement" : "complement"
},
   "phone" : "phone",
   "name" : "name",
   "id" : 0,
   "email" : "email",
   "dni" : "dni"
}
```

### **Produces**

This API call produces the following media types according to the Accept request header; the media type will be conveyed by the Content-Type response header.

• \*/\*

# Responses

200

**OK DataPatient** 

DELETE /api/patient/{id}



(deactivate Patient)

### Path parameters

```
id (required)
```

Path Parameter - format: int64

## Return type

String

### Example data

Content-Type: application/json

```
пп
```

### **Produces**

This API call produces the following media types according to the Accept request header; the media type will be conveyed by the Content-Type response header.

• \*/\*

# Responses

200

**OK String** 

```
GET /api/patient
(getAll)
```

# <u>Up</u>

### **Query parameters**

pagination (required) Query Parameter —

## Return type PageDataPatientList

# Example data

Content-Type: application/json

```
"number" : 5,
"size" : 1,
"last" : true,
"numberOfElements" : 2,
"totalPages" : 0,
"pageable" : {
  "paged" : true,
  "pageNumber" : 9,
  "offset": 7,
  "pageSize" : 3,
  "unpaged" : true
"sort" : {
  "unsorted" : true,
  "sorted" : true,
  "empty" : true
"content" : [ {
  "name" : "name",
  "id" : 5,
  "email" : "email",
  "dni" : "dni"
}, {
  "name" : "name",
```

```
"id" : 5,
    "email" : "email",
    "dni" : "dni"
} ],
    "first" : true,
    "totalElements" : 6,
    "empty" : true
}
```

### **Produces**

This API call produces the following media types according to the Accept request header; the media type will be conveyed by the Content-Type response header.

• \*/\*

# Responses

200

**OK PageDataPatientList** 

```
GET /api/patient/{id}
(getOne)

Path parameters
  id (required)
    Path Parameter — format: int64
```

# Return type Patient

# Example data

Content-Type: application/json

```
{
   "address" : {
        "zipcode" : "zipcode",
        "address" : "address",
        "city" : "city",
        "county" : "county",
        "complement" : "complement"
   },
   "phone" : "phone",
   "name" : "name",
   "id" : 6,
   "isActive" : true,
   "email" : "email",
   "dni" : "dni"
}
```

#### **Produces**

This API call produces the following media types according to the Accept request header; the media type will be conveyed by the Content-Type response header.

• \*/\*

# Responses

200

**OK** Patient

<u>Up</u>

# POST /api/patient/search



### (searchByCriteria)

### Consumes

This API call consumes the following media types via the Content-Type request header:

• application/json

### Request body

```
body patient_search_body (required)
Body Parameter —
```

# Return type

**PageDataPatientList** 

### Example data

Content-Type: application/json

```
"number" : 5,
  "size" : 1,
  "last" : true,
  "numberOfElements" : 2,
  "totalPages" : 0,
  "pageable" : {
    "paged" : true,
    "pageNumber" : 9,
    "offset": 7,
    "pageSize" : 3,
    "unpaged" : true
  "sort" : {
    "unsorted" : true,
    "sorted" : true,
    "empty" : true
  "content" : [ {
    "name" : "name",
    "id" : 5,
    "email" : "email",
    "dni" : "dni"
  }, {
    "name" : "name",
    "id" : 5,
    "email" : "email",
    "dni" : "dni"
  } ],
  "first" : true,
  "totalElements" : 6,
  "empty" : true
}
```

### **Produces**

This API call produces the following media types according to the Accept request header; the media type will be conveyed by the Content-Type response header.

• \*/\*

# Responses

200

**OK PageDataPatientList** 

PUT /api/patient



(updateDoctor)

#### Consumes

This API call consumes the following media types via the Content-Type request header:

• application/json

## Request body

body <u>DataPatientUpdate</u> (required) Body Parameter —

### Return type

String

### Example data

Content-Type: application/json

...

### **Produces**

This API call produces the following media types according to the Accept request header; the media type will be conveyed by the Content-Type response header.

• \*/\*

### Responses

200

OK **String** 

# UserAuthentication

POST /api/auth/login



(authenticateUser)

#### Consumes

This API call consumes the following media types via the Content-Type request header:

• application/json

### Request body

body **DataUserAuthentication** (required)

Body Parameter –

### Return type

String

### Example data

Content-Type: application/json

```
пп
```

### **Produces**

This API call produces the following media types according to the Accept request header; the media type will be conveyed by the Content-Type response header.

• \*/\*

# Responses 200

**OK String** 

POST /api/auth/register

<u>Up</u>

(createUser)

### Consumes

This API call consumes the following media types via the Content-Type request header:

• application/json

### Request body

```
body <u>DataUserRegister</u> (required)
Body Parameter —
```

### Return type DataUser

### Example data

Content-Type: application/json

```
{
   "password" : "password",
   "id" : 0,
   "username" : "username"
}
```

### **Produces**

This API call produces the following media types according to the Accept request header; the media type will be conveyed by the Content-Type response header.

• \*/\*

### Responses

200

**OK DataUser** 

# Models

[ Jump to Methods ]

### **Table of Contents**

1. Address

2	<u>Appointment</u>	
	<u>DataAddress</u>	
	<u>DataAppointment</u>	
	<u>DataCancelAppointment</u>	
	<u>DataDoctor</u>	
	<u>DataDoctorList</u>	
	<u>DataDoctorRegister</u>	
9.	<u>DataDoctorUpdate</u>	
	<u>DataPatient</u>	
	<u>DataPatientList</u>	
	<u>DataPatientRegister</u>	
	<u>DataPatientUpdate</u>	
	<u>DataScheduleAppointment</u>	
	DataUser  DataUser	
	<u>DataUserAuthentication</u> DataUserBegisten	
	<u>DataUserRegister</u> <u>Doctor</u>	
	<u>PageDataDoctorList</u>	
	<u>PageDataPatientList</u>	
	Pageable	
	PageableObject	
	<u>Patient</u>	
	<u>SortObject</u>	
25.	<pre>doctor_search_body</pre>	
26.	<pre>patient_search_body</pre>	
Ad	dress	<u>Up</u>
	address (optional)	
	<u>String</u>	
	county (optional)	
	<u>String</u>	
	-	
	zipcode (optional) <u>String</u>	
	•	
	complement (optional)	
	<u>String</u>	
	city (optional)	
	<u>String</u>	
۸n	pointment	<u>Up</u>
ΑÞ	pointment	<u> </u>
	id (optional)	
	<u>Long</u> format: int64	
	doctor (optional)	
	<u>Doctor</u>	
	patient (optional)	
	<u>Patient</u>	
	date (optional)	
	<u>Date</u> format: date-time	
	cancellationReason (optional)	
	<u>String</u>	
	Enum:	
	PATIENT_CANCELED DOCTOR_CANCELED	
	DOCTON_CANCELED	

OTHER

DataAddress	
address <u>String</u>	
county <u>String</u>	
zipcode <u>String</u>	
complement <u>String</u>	
city <u>String</u>	
DataAppointment	<u>Up</u>
id (optional) <u>Long</u> format: int64	
patient (optional) <u>Patient</u>	
doctor (optional) <u>Doctor</u>	
date (optional) <u>Date</u> format: date-time	
DataCancelAppointment	<u>Up</u>
idAppointment <u>Long</u> format: int64	
cancellationReason (optional) <u>String</u>	
Enum:  PATIENT_CANCELED  DOCTOR_CANCELED  OTHER	
DataDoctor	<u>Up</u>
id (optional) <u>Long</u> format: int64	
name (optional) <u>String</u>	
email (optional) <u>String</u>	
phone (optional) <u>String</u>	
dni (optional) String	

address (optional) <u>DataAddress</u> DataDoctorList <u>Up</u> id (optional) Long format: int64 name (optional) <u>String</u> email (optional) <u>String</u> dni (optional) **String** specialty (optional) <u>String</u> Enum: **ORTHOPEDY CARDIOLOGY PEDIATRY DERMATOLOGY ONCOLOGY** DataDoctorRegister <u>Up</u> name <u>String</u> email <u>String</u> phone String dni **String** specialty <u>String</u> Enum: **ORTHOPEDY CARDIOLOGY PEDIATRY DERMATOLOGY** ONCOLOGY address **DataAddress** DataDoctorUpdate <u>Up</u> id Long format: int64 name (optional) <u>String</u> dni (optional)

**String** 

address (optional) <u>DataAddress</u>

DataPatient	
id (optional) <u>Long</u> format: int64	
name (optional) <u>String</u>	
email (optional) <u>String</u>	
phone (optional) <u>String</u>	
dni (optional) <u>String</u>	
address (optional) <u>DataAddress</u>	
DataPatientList	<u>Up</u>
id (optional) <u>Long</u> format: int64	
name (optional) <u>String</u>	
email (optional) <u>String</u>	
dni (optional) <u>String</u>	
DataPatientRegister	
name <u>String</u>	
email <u>String</u>	
phone <u>String</u>	
dni <u>String</u>	
address <u>DataAddress</u>	
DataPatientUpdate	
id <u>Long</u> format: int64	
name (optional) <u>String</u>	

dni (optional) <u>String</u> address (optional) **DataAddress DataScheduleAppointment** <u>Up</u> idPatient Long format: int64 idDoctor (optional) Long format: int64 date **Date** format: date-time specialty (optional) <u>String</u> Enum: **ORTHOPEDY CARDIOLOGY PEDIATRY DERMATOLOGY** ONCOLOGY DataUser <u>Up</u> id (optional) Long format: int64 username (optional) <u>String</u> password (optional) <u>String</u> <u>Up</u> **DataUserAuthentication** username (optional) <u>String</u> password (optional) <u>String</u> <u>Up</u> DataUserRegister username <u>String</u> password <u>String</u> <u>Up</u> **Doctor** id (optional) Long format: int64

	name (optional) <u>String</u>	
	email (optional) <u>String</u>	
	phone (optional) String	
	dni (optional) <u>String</u>	
	isActive (optional) <u>Boolean</u>	
	specialty (optional) <u>String</u> Enum:  ORTHOPEDY  CARDIOLOGY  PEDIATRY  DERMATOLOGY  ONCOLOGY	
	address (optional) <u>Address</u>	
Pā	ngeDataDoctorList	<u>Up</u>
	totalPages (optional) <a href="mailto:lnteger">lnteger</a> format: int32	
	totalElements (optional) <u>Long</u> format: int64	
	size (optional) <a href="mailto:lnteger">lnteger</a> format: int32	
	content (optional) <pre>array[DataDoctorList]</pre>	
	number (optional) <a href="mailto:lnteger">lnteger</a> format: int32	
	sort (optional) <u>SortObject</u>	
	first (optional) <u>Boolean</u>	
	last (optional) <u>Boolean</u>	
	numberOfElements (optional) <a href="mailto:lnteger">lnteger</a> format: int32	
	pageable (optional) <u>PageableObject</u>	
	empty (optional) <u>Boolean</u>	
Pa	ngeDataPatientList	<u>Up</u>
	totalPages (optional) <a href="mailto:lnteger">lnteger</a> format: int32	

 $file: ///C: /Users/juanj/AppData/Local/Temp/Temp1\_html-client-generated.zip/index.html$ 

totalElements (optional) **Long** format: int64 size (optional) *Integer* format: int32 content (optional) array[DataPatientList] number (optional) *Integer* format: int32 sort (optional) <u>SortObject</u> first (optional) **Boolean** last (optional) **Boolean** numberOfElements (optional) *Integer* format: int32 pageable (optional) <u>PageableObject</u> empty (optional) **Boolean** Pageable <u>Up</u> page (optional) <u>Integer</u> format: int32 size (optional) **Integer** format: int32 sort (optional) array[String] PageableObject <u>Up</u> offset (optional) **Long** format: int64 sort (optional) <u>SortObject</u> pageNumber (optional) <u>Integer</u> format: int32 pageSize (optional) **Integer** format: int32 unpaged (optional) **Boolean** paged (optional) Boolean **Patient** <u>Up</u> id (optional)

Long format: int64 name (optional) <u>String</u> email (optional) <u>String</u> phone (optional) <u>String</u> dni (optional) <u>String</u> isActive (optional) <u>Boolean</u> address (optional) **Address** SortObject <u>Up</u> empty (optional) <u>Boolean</u> sorted (optional) <u>Boolean</u> unsorted (optional) **Boolean** doctor\_search\_body <u>Up</u> dataDoctorList (optional) <u>DataDoctorList</u> pagination (optional) <u>Pageable</u> patient\_search\_body <u>Up</u> dataPatientList (optional) **DataPatientList** pagination (optional) <u>Pageable</u>